

HDL-CHOLESTEROL

| CONTENTS | | | |
|---|---------|-----------------|-----------|
| REF | 1133010 | HDL-Cholesterol | 2 x 40 mL |
| For <i>in vitro</i> diagnostic use only | | | |

HDL-CHOLESTEROL

DIFFERENTIAL PRECIPITATION

Enzymatic colorimetric test

ENDPOINT

PRINCIPLE

This technique¹ uses a separation method based on the selective precipitation of apolipoprotein B-containing lipoproteins (VLDL, LDL and (a)Lpa) by phosphotungstic acid/MgCl₂, sedimentation of the precipitant by centrifugation, and subsequent enzymatic analysis of high density lipoproteins (HDL) as residual cholesterol remaining in the clear supernatant.


REAGENT COMPOSITION

R1 **Precipitating reagent.** Phosphotungstic acid 0.63 mmol/L, magnesium chloride 25 mmol/L. Stabilizers.

CAL **HDL-Cholesterol standard.** Cholesterol 50 mg/dL (0.38 mmol/L). Organic matrix based primary standard.

Cholesterol MR. Optative. Ref: 1118005, 1118010, 1118015.

STORAGE AND STABILITY

 Store at 2-8°C.
The reagents are stable until the expiry date stated on the label.

REAGENT PREPARATION

The Reagents and Standard are ready-to-use.

SAMPLES

Serum or EDTA plasma free of hemolysis obtained by the patient after an overnight fast. Remove from cells within 3 hours of venipuncture. Samples may be kept at 4-8°C for 2 weeks, and at -20°C for 3 months with no alteration of HDL cholesterol. The supernate containing the HDL fraction is conveniently prepared on the day of sample collection and may be analysed after 2 weeks at 4-8°C or 3 months at -20°C in a non-selfdefrosting freezer.²

INTERFERENCES

- Hemoglobin (>200 mg/dL) and bilirubin (>10 mg/dL) do not interfere with the test.³
- Turbidity in samples may indicate elevated triglycerides or non-fasting specimens.

MATERIALS REQUIRED

I. Precipitation

- Dilutor and pipettes.
- Centrifuge tubes (13 x 100 m/m).
- Vortex mixer.
- Desktop centrifuge.

II. Colorimetry

- Kit for measuring Total Cholesterol.
- Constant temperature incubator set at 37°C.
- Photometer or colorimeter capable of measuring absorbance at 550 ± 10 nm.

PROCEDURE

I. Precipitation

1. Bring reagents and samples to room temperature.
2. Pipette into labelled centrifuge tubes:

| | | |
|-----------------------|--------|--|
| Sample or Standard | 0.2 mL | Ratio $\frac{\text{Sample}}{\text{Reagent}} = \frac{1}{2}$ |
| Precipitating reagent | 0.4 mL | |
| | | Dil. factor = 3 |

3. Vortex and allow to stand for 10 minutes at room temperature.
4. Centrifuge for 10 minutes at 4000 r.p.m., or 2 minutes at 12000 r.p.m.
5. Separate off the clear supernatant within 2 hours.

In case of turbid supernatants caused by elevated triglycerides (>350 g/dL) the sample should be diluted 1:2 with saline and steps 2,3,4 and 5 repeated. Multiply the result of the colorimetry by 2.

II. Colorimetry

1. Bring the Cholesterol MR Monoreagent and the cholesterol standard (50 mg/dL) of the kit to room temperature
2. Pipette into labelled tubes:

| TUBES | Blank | Sample Supernat | Standard Supernat |
|-------------|--------|-----------------|-------------------|
| Monoreagent | 1.0 mL | 1.0 mL | 1.0 mL |
| Supernat | – | 50 µL | – |
| Standard | – | – | 50 µL |

3. Mix and let stand the tubes for 10 minutes at room temperature or 5 minutes at 37°C.
4. Read the absorbance (A) of the supernatant and the standard at 550 nm against the reagent blank.

The color is stable for at least 30 minutes protected from light.

CALCULATIONS

$$\frac{A_{\text{Supernatant}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = \text{mg/dL HDL-Cholesterol}$$

If results are to be expressed as SI units apply:
 $\text{mg/dL} \times 0.0259 = \text{mmol/L}$

REFERENCE VALUES⁴

Clinical values of HDL-Cholesterol used to classify risk groups.

| Cholesterol from lipoproteins of high density | | RISK |
|---|--------------------------------|----------|
| Men | > 55 mg/dL (> 1.42 mmol/L) | Low |
| | 35-55 mg/dL (0.90-1.42 mmol/L) | Moderate |
| | < 40 mg/dL (< 1.04 mmol/L) | High |
| Women | > 65 mg/dL (> 1.68 mmol/L) | Low |
| | 45-65 mg/dL (1.16-1.68 mmol/L) | Moderate |
| | < 45 mg/dL (< 1.16 mmol/L) | High |

CLINICAL SIGNIFICANCE

Low HDL-cholesterol is a strong independent predictor of coronary heart disease. In ATP III⁴, low HDL cholesterol is defined categorically as a level < 40 mg/dL (1.04 mmol/L), a change from the level of < 35 mg/dL in ATP II (1993).

Low HDL cholesterol is used as a risk factor to estimate 10-year risk for coronary heart disease, having several causes: elevated triglycerides, overweight and obesity, physical inactivity, and type 2 diabetes. Other causes are, cigarette smoking, very high carbohydrate intakes (> 60% of calories), and certain drugs as anabolic steroids and progestational agents.

ANALYTICAL PERFORMANCE

- **Linearity.** Up to 275 mg/dL

- **Precision**

| mg/dL | Within-run | | | Between-run | | |
|-------|------------|------|------|-------------|------|------|
| Mean | 42.1 | 45.8 | 54.6 | 42.1 | 45.8 | 54.6 |
| SD | 0.23 | 0.23 | 0.2 | 0.27 | 0.28 | 0.31 |
| CV% | 0.54 | 0.5 | 0.34 | 0.64 | 0.61 | 0.52 |
| N | 10 | 10 | 10 | 10 | 10 | 10 |

Replicates: 10 for each level.

Replicates: 10 for each level

Instrument: CECIL CE 2021

for 8 days.

- **Sensitivity.** Using a 1:3 sample/reagent at 550 nm, 10 mg of cholesterol will produce a net absorbance of approximately 0.037.

- **Correlation.** This assay (y) was compared with a similar commercial method (x). The results were:

$$N = 25 \quad r = 0.995 \quad y = 0.985 + 2.6$$

REFERENCES

1. Burstein, M., Scholnick, H.R. and Morfin, R. Scand. J. Clin. Lab. Invest. 40 : 560 (1980).
2. Finley, P.R., Shifman, R.B., Williams, R.S. and Lichti, D.I. Clin. Chem. 24 : 931 (1971).
3. Tietz, N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).