MULTICALIBRATOR CC/H 

PRESENTACION / CONTENTS

REF	1975005	Multicalibrator CC/H	5 x 5 mL
-----	---------	----------------------	----------

Sólo para uso diagnóstico *in vitro*
Only for *in vitro* diagnostic use

MULTICALIBRADOR CC/H
MULTICALIBRATOR CC/H

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INTENDED USE

LINEAR Multicalibrator CC/H is a lyophilized human serum for use in clinical chemistry assays. This calibrator is offered in the form of an assayed serum. Given values are assigned on results from 3000 different laboratories. The concentrations and activities have been selected so as to ensure optimum calibration for use in manual and automatic analyzers.

COMPOSITION

CAL Prepared from human serum, no reactive conservatives and stabilizers.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV 1/2. However handle cautiously as potentially infectious.

STABILITY

Unreconstituted serum is stable up to the expiry date shown on the label when stored at 2-8°C. Once reconstituted, are stable for at least 8 hours at 25°C, 7 days at 2-8°C or 1 month at -20°C when frozen once (see Notes).

RECONSTITUTION

1. Open the vial carefully, avoiding any loss of material.
2. Reconstitute by pipetting exactly **5 mL** of distilled water at 20-25°C, into the vial.
3. Replace the rubber stopper and leave to stand for 30 minutes out of bright light before using.
4. Swirl gently several times during the reconstitution period to ensure that the contents are completely dissolved.
5. Prior to use, mix the contents by inverting the vial. Do not shake the vial as the formation of foam should be avoided.
6. The **Multicalibrator CC/H** is then ready for use.

NOTES

- The bilirubin present in the serum is light sensitive and it is recommended that the serum be stored in the dark. Stored in the dark it is stable for at least 8 hours at 2-8°C. Do not store at 15-25°C. Do not freeze. For Total and Prostatic acid phosphatase the material should be stabilized by adding 1 drop (25-30 µL) of 0.7M acetic acid solution to 1mL of serum. After stabilization Total and Prostatic acid phosphatase are stable for at least 2 hours at 25°C, at least 2 days at 2-8°C and 1month when frozen once at -20°C. Alkaline phosphatase levels in the reconstituted serum will rise over the stability period. It is recommended that the reconstituted serum be allowed to stand for 1 hour at 25°C before measurement.
- Bacterial contamination of the **Multicalibrator CC/H** will cause reductions in the stability of many components.
- Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.
- The assigned values are lot-specific. The exact values are given in the value sheet.

MULTICALIBRATOR CC/H

REF 1975005		LOT N° 19758		EXP 2018-12
COMPONENTE COMPONENT	VALOR* VALUE*	UNIDAD UNIT	MÉTODO METHOD	
Acid Phosphatase (non-prostatic)	14.77	U/L	1-Naphthyl Phosphate, Kinetic with Pentane diol Activation 37°C	
Acid Phosphatase (prostatic)	39.23	U/L	1-Naphthyl Phosphate, Kinetic with Pentane diol Activation 37°C	
Acid Phosphatase (total)	54.0	U/L	1-Naphthyl Phosphate, Kinetic with Pentane diol Activation 37°C	
Albúmina Albumin	29.4 2.94	g/L g/dL	Bromocresol Green (BCG)	
Fosfatasa Alcalina Alkaline Phosphatase	456 356 292	U/L	Diethanolamine buffer DEA 37°C Diethanolamine buffer DEA 30°C Diethanolamine buffer DEA 25°C	
ALT (GPT)	120 89 68	U/L	Tris buffer no P5P IFCC 37°C Tris buffer no P5P IFCC 30°C Tris buffer no P5P IFCC 25°C	
α -Amilasa α -Amylase	282	U/L	2-chloro-pNPG3 37°C	
AST (GOT)	142 97 68	U/L	Tris buffer no P5P IFCC 37°C Tris buffer no P5P IFCC 30°C Tris buffer no P5P IFCC 25°C	
Bilirrubina Directa Bilirubin Direct	29.5 1.73	μ mol/L mg/dL	Modified Jendrassik & Grof. DIAZOTIZED SULFANILIC	
Bilirrubina Total Bilirubin Total	86.1 5.04	μ mol/L mg/dL	Modified Jendrassik & Grof. DIAZOTIZED SULFANILIC	
Calcio Calcium	3.1 12.3	mmol/L mg/dL	Arsenazo III	
Calcio Calcium	3.1 12.4	mmol/L mg/dL	Cresolphthalein complexone	
Cloruros Chloride	114	mmol/L	THIOCYANATE	
Colesterol Cholesterol	7.15 276	mmol/L mg/dL	CHOLESTEROL OXIDASE/PEROXIDASE	
Colinesterasa Cholinesterase	5213	U/L	Butyrylthiocholine 37°C	
CK Total	508 318 217	U/L	CK NAC-activated (IFCC) 37°C CK NAC-activated (IFCC) 30°C CK NAC-activated (IFCC) 25°C	
Creatinina Creatinine	351 3.97	μ mol/L mg/dl	Alkaline picrate with deproteinization / Creatinine E Enzymatic	

*Mean of instruments / Media de los instrumentos

MULTICALIBRATOR CC/H

REF 1975005	LOT N° 19758	EXP 2018-12
--------------------	---------------------	--------------------

COMPONENTE COMPONENT	VALOR* VALUE*	UNIDAD UNIT	MÉTODO METHOD
Glucosa Glucose	15.21 274	mmol/L mg/dL	Glucose oxidase/ Peroxidase
γ-GT	136 107 85	U/L	Gamma glutamyl.-3-carboxy-4-nitroanilide 37°C Gamma glutamyl.-3-carboxy-4-nitroanilide 30°C Gamma glutamyl.-3-carboxy-4-nitroanilide 25°C
Hierro Iron	34.9 195	μmol/L μg/dL	Ferrozine.
Hierro Iron	35.44 198	μmol/L μg/dL	Cromazurol S
LD (LDH)	639 462 324	U/L	P->L SFBC 37°C P->L SFBC 30°C P->L SFBC 25°C
Magnesio Magnesium	1.68 4.08	mmol/L mg/dL	Calmagite
Fósforo Phosphorus	2.23 6.91	mmol/L mg/dL	Phosphomolybdate UV
Potasio Potassium	5.86	mmol/L	Pyruvate kinase
Proteínas totales Total protein	55 5.50	g/L g/dL	Biuret reaction end point
Sodio Sodium	158	mmol/L	Sodium-dependent β-galactosidase
TIBC	47.4 265	μmol/L μg/dL	FE+UIBC (saturation with iron)
Triglicéridos Triglycerides	2.78 246	mmol/L mg/dL	Lipase/Glycerol Dehydrogenase
Urea	17.4 105	mmol/L mg/dL	Berthelot Urease/Salicylate
Urea	18.0 108	mmol/L mg/dL	BUN Urease/GIDH
Ácido Úrico Uric Acid	572 9.61	μmol/L mg/dL	Uricase / peroxidase

*Mean of instruments / Media de los instrumentos

B1975-2/0509
R1