

# qLabs<sup>®</sup> ElectroMeter Plus

for Professional Point of Care



MICROPOINT

# qLabs® ElectroMeter Plus

for Professional Point of Care

## Diagnostics on Demand for POC

### Handheld Rapid Testing Platform

#### Features

- On board quality controls
- Fingertip blood sample of 10 µL
- Reliable laboratory quality results in less than 7 minutes
- External liquid QC (LQC) available
- High isolation between channels for high precision
- Built-in temperature controlled testing
- eStation for printing and connectivity
- Built-in barcode scanner
- No calibration or instrument manipulation needed

#### Current Test Menu

- PT
- aPTT
- PT/aPTT Combo

#### ABOUT MICROPOINT

Micropoint Bioscience, Inc. is a global provider of Point of Care Testing (POCT) systems and solutions. We are improving patient care by providing rapid, accurate, reliable, and low cost diagnostic products. Micropoint's patented microfluidic technologies and three proprietary POCT platforms enable testing at the same high quality level as a central lab, but with improved speed and convenience. Micropoint's products are manufactured in cGMP and ISO 13485 certified facilities.

<b>Interface</b> Display Barcode scanner	2.1 inch 128x96 LCD Laser class 2 according to EN 60825-1:2007
<b>Software</b> Language	Multi-language menu
<b>Memory Capacity</b> Patient records	Stores up to 200 patient records and 12 QC results
<b>Performance</b> Measurable Range	INR between 0.5 to 7.5 APTT from 20 sec to 420 sec
<b>Specifications</b> Dimensions Weight Operating conditions  Power supply	135mm x 65mm x 34.5mm (L x W x H) 154 g (Without Batteries) Temperature: 10 °C - 35 °C Relative humidity: 10 % - 90 % 100V-240V AC 4 X AA batteries



Micropoint Bioscience, Inc.  
 3333 Octavius Drive, Suite 110  
 Santa Clara, CA 95054  
 Tel : +1-408-588-1682  
 Fax : +1-408-588-1620  
 E-mail : info@micropointbio.com  
 www.micropointbio.com